


Airmonte A5000 Alternating Air Overlay



Airmonte A5000 Alternating Air Overlay


The Airmonte A5000 is an intuitive advanced alternating air overlay. Trusted within care facilities across Australia, the Airmonte A5000 delivers a high standard and consistency of care, with a flawless, easy to operate advanced design at excellent value. Key clinical requirements including adjustable cell pressure, alternating/static modes, adjustable comfort level, quick-release CPR valves, quality componentry, effortless operation and easy maintenance ensures effectiveness and comfort for patients and support for carers. Select the Airmonte Pressure Care Mattresses and achieve unsurpassed reliability and exceptional care.



/ Airmonte A5000

CODE	SIZE (MM)	
A5000	1960x860x125	Single
A5KS1000		King Single
THERAPEUTIC SUPPORT LOADING*	25 - 150KG	
PRODUCT WEIGHT	10KG	
WARRANTY**	12 MONTH WARRANTY	
APPLICATIONS	✓ AGED CARE ✓ HOME CARE	

Features



ALTERNATING CELL CONSTRUCTION

- ✓ 24 Interchangeable alternating cells providing excellent pressure re-distribution
- ✓ Patient is supported constantly, simultaneously allowing oxygen flow to all critical areas, preventing and/or healing pressure sores
- ✓ Height of mattress (125mm) allows maximum patient immersion for effectiveness and comfort

MODULAR DESIGN

- ✓ Every component is modular and replaceable for low ongoing operational costs

UNIQUE QUICK-CONNECT SYSTEM

- ✓ Enables simple pump set-up
- ✓ Quick and easy for caregiver operation

QUALITY, ERGONOMIC PUMP UNIT

- ✓ Styling is compact and efficient, ergonomically designed with an easy to operate panel
- ✓ Quality componentry

PUMP CONTROL PANEL FUNCTIONS

- ✓ Silent operation minimizing patient disturbance
- ✓ Intuitive, simple operation
- ✓ One-touch quick-change air filter for ease of maintenance

EASY TO ACCESS CPR


- ✓ Positioned on air supply tubes

FIRM FOAM SURROUND

- ✓ for mattress dimensions greater than 1960x860x125mm
- ✓ Firm, Premium Grade Sides
- ✓ Supportive Foundation Foam: 25mm Thick

Airmonte A5000 Pressure Mapping

SUPINE



NOTES

Pressure Mapping Parameters

Bed Surface: Alrick, 5 Function Electric Hospital/Homecare Bed

Ambient Temperature: 21°C

Duration: Pressure Mapping was captured after 6 hours of Subject Laying Upon Support Surface

Subject: 98kg, Male, Height: 188cm



P. 1300 73 11 63

W. fortehealthcare.com.au

E. sales@fortehealthcare.com.au

