




Custom Printing



MEDICAL MATTRESS EXPERTS

Forté Product Printing & Audit Recording Customisation

All Forté Healthcare products are printed with an indelible digital print which meets all TGA requirements and includes other relevant information for the user.

Forté Healthcare has an in-house design and printing service and can customise printing including custom branding or facility specific information such as ward name, as an example.

TGA PRINT REQUIREMENTS

- Product Dimensions
- Product Weight
- Product Name/Description
- Therapeutic Weight Range
- Cleaning & Care Instructions
- Manufacturers Details
- Date of Manufacture

COMMON PRINT INFORMATION FOR USERS

- Heel slope
- Mattress positioning direction
- Carry handle location
- Evacuation system functions

COMMON CUSTOM PRINTS

- Product tracking and control information i.e. Ward Locations, Facility Name, Asset Number
- Additional Cleaning & Care information specific to your facility
- Custom branding
- Audit Systems/records specific to your facility

AVAILABLE ON ANY FORTÉ MEDICAL COVER





P. 1300 73 11 63

W. fortehealthcare.com.au

E. sales@fortehealthcare.com.au