




Forté Product Customisations



MEDICAL MATTRESS EXPERTS

Forté Firm Safety Sides Customisation

A Firm Safety Perimeter is an accessory that can be added to any Forté mattress to improve patient safety by reducing risk of entrapment and providing enhanced mattress stability – especially pediatric patients. This unobtrusive firm perimeter is made from a high-density foam core and has in-built hinges on the sides to allow for the mattress to easily conform to varied bed positions and without causing foam fatigue. Standard Thickness: 75mm.

SUBTLE PATIENT RESTRAINT


- Manufactured with quality foam which conforms to patient transfer and primary care whilst maintaining its profile.
- Unobtrusive and does not impede primary patient care or general bedside nursing duties.
- Comfortable and protective
- Provides a non-restraining and secure feeling
- Encourages independent repositioning and exits from the bed
- Prevents risk of entrapment and falls

ENCOURAGES CENTRAL PATIENT POSITIONING

- Prevents patients rolling or falling off the mattress and reduces risk of entrapment for restless and agitated or pediatric patients
- Central patient positioning ensures full therapeutic function of the mattress

AVAILABLE WITH ANY FORTÉ MATTRESS AND COVER

- Must be purchased in conjunction with a Forté mattress
- [Read more about Our range of mattresses](#)




Concave Safety Sides

Concave Safety Sides are an accessory that can be added to any Forté mattress for those with falls risk. These raised sides encourage central patient positioning and gentle protection for patients included to turn out of bed.

Forté Healthcare's Concave Safety Sides are manufactured with quality foam which confirms to patient transfer and primary care whilst maintaining its profile.

Standard Thickness: 95mm, Standard Height/Slope: 245-150mm. Customisable

- ✓ **FALLS RISK PREVENTION**
- ✓ **ENCOURAGES CENTRAL PATIENT POSITIONING**
- ✓ **SUBTLE PATIENT RESTRAINT**



Forté Mattress Carry Handles

Carry Handles are an accessory that can be added to the side of any Forté mattress to assist with mattress handling.

Forté Healthcare's unique Wipe-Clean, non-absorbant Carry Handles meet all infection control concerns. The material used is sealed and all attachment points are high pressure fastened (not sewn) preventing ingress and penetration of moisture and fluid. Forté Carry Handles can be wiped cleaned. Often in the course of use, infectious fluid will spill down the side of the mattress cover but with this handle, nothing can soak in and contaminate the carry point.

The handle system also serves to assist with mattress handling and prevents damage to mattress cover caused by general mishandling.

A unique Flush mount design averts catching on ancillary/nearby patient care equipment and elements of the bed frame, decreasing the risk of damage to the mattress cover

- ✓ **MEETS WH&S AND INFECTION CONTROL CONCERNS**
- ✓ **STURDY**



Forté Product Printing & Audit Recording Customisation

All Forté Healthcare products are printed with an indelible digital print which meets all TGA requirements and includes other relevant information for the user.

Forté Healthcare has an in-house design and printing service and can customise printing including custom branding or facility specific information such as ward name, as an example.

TGA PRINT REQUIREMENTS

- Product Dimensions
- Product Weight
- Product Name/Description
- Therapeutic Weight Range
- Cleaning & Care Instructions
- Manufacturers Details
- Date of Manufacture

COMMON PRINT INFORMATION FOR USERS

- Heel slope
- Mattress positioning direction
- Carry handle location
- Evacuation system functions

COMMON CUSTOM PRINTS

- Product tracking and control information i.e. Ward Locations, Facility Name, Asset Number
- Additional Cleaning & Care information specific to your facility
- Custom branding
- Audit Systems/records specific to your facility

AVAILABLE ON ANY FORTÉ MEDICAL COVER



Forté Product Sizing & Dimesion Customisation

Every patient is unique and no two care environments are necessarily the same.

We readily acknowledge that with patients well-being and staff safety is at stake, we cannot take a 'one-size suits all' approach.

All Forté Healthcare products are manufactured and designed in the state of the art manufacturing facility in Armidale, Australia. Common sizes and models are held in stock, the rest is made to order. This means Forté can customize any product to suit your specific needs.

Sizing & Dimension Customizations

Any standard or custom product can be made to suit any requirement and bed platform – from neonatal and paediatric needs through to king size luxury domestic beds, truncated stretcher mattress or complex acute care bed systems. No longer do caregivers and patients have to compromise on achieving the perfect fit and function and we can ensure comfort and aesthetic requirements are met. There is so much depending upon the Sleep Surrface performing correctly and our bespoke, individualized approach ensures optimal outcomes.

Some of the favourite size customisations include:

Split Queens

Double/Queen/King Size Zephair Air Mattress System Partner mattress

Cuddle Beds (read more below)

Extra Height - 250-350mm increased height to match existing bed frames or achieve optimal transfer height

Rest assured, All customisations are checked off by our Research & Development team to ensure the therapeutic function of the mattress is not impeded by the requested customisation.

Speak to Forté today to request further information!

Cuddle Bed Mattresses

Whilst Forté Healthcare designs mattresses specifically for pressure care, we have a strong emphasis on patient well-being, comfort, compliance and holistic care.

Much research shows that offering a high level of care whilst accommodating as much independence and normality as possible has powerful, positive impacts on people's lives.

Forté will customise any mattress to accompany a 'Cuddle Bed', encouraging normal partner sleeping and minimising anxiety, isolation, strangeness among patients. This includes in acute care applications such as palliative and maternity wards.

Forté Healthcare offer 'Cuddle Bed' solutions in low-high risk static foam surfaces and also high-very high risk air systems

Speak to Forté today ask what can be done to help!



P. 1300 73 11 63

W. fortehealthcare.com.au

E. sales@fortehealthcare.com.au