




SONNET PRESSURE CARE POSITIONING DEVICES

Sonnet Heel Elevator



The Forté Heel Elevator supports the patients leg and providing suspension to the critical heel region. Manufactured using premium immersion foams and a pressure care covering, the Forté Heel Wedge provides effective support for a wide range of patients.



Standard Cover:

ADVANTIFLEX

PROFILED PRESSURE RELIEVING CORE

- ✓ High quality, open cell foam
- ✓ Anti-microbial treatment
- ✓ Tapered for maximum patient comfort

AIROFORM MEMORY IMMERSION INTERFACE LAYER

- ✓ Pressure sensitive
- ✓ High temperature stability
- ✓ Visco-elastic cell structure provides excellent immersion and envelopment
- ✓ Comfort: Soft

/ Sonnet Heel Elevator

CODE	SIZE (MM)
HE1000	400mm long

Sonnet Heel Cushion



The Forté Pressure Relieving Heel Cushion Provides effective protection to the critical Heel region. Manufactured to be able to be used in a variety positions, the Forté Heel Cushion requires minimal education for clinical care staff and is popular within acute and frail aged care for its simplistic design and effectiveness. Manufactured using premium immersion foam and a pressure care covering, the Forté Heel Cushion is suited for a wide range of patients.



/ Sonnet Heel Cushion

CODE	SIZE (MM)
HC1000	670mm wide

Standard Cover:
ADVANTIFLEX


- PRESSURE RELIEVING CORE**
- ✓ High quality, open cell foam
 - ✓ Anti-microbial treatment
 - ✓ Simplistic 4-way turn positioning

- AIROFORM MEMORY IMMERSION INTERFACE**
- ✓ Pressure sensitive
 - ✓ High temperature stability
 - ✓ Visco-elastic cell structure provides excellent immersion and envelopment

Sonnet Heel Wedge



The Sonnet Heel Wedge is used for leg elevation and heel suspension for effective pressure relief to critical heel areas. Manufactured using premium immersion foams and a pressure care covering, the Forté Heel Wedge is suited for a wide range of patients and ideal for performance and longevity.



/ Sonnet Heel Wedge


CODE	SIZE (MM)
HW1000	280mm wide
HW1001	400mm wide
HW1002	680mm wide

- Standard Cover:
ADVANTIFLEX
- AIROFORM MEMORY IMMERSION INTERFACE**
- ✓ Pressure sensitive
 - ✓ High temperature stability
 - ✓ Visco-elastic cell structure provides excellent immersion and envelopment
- SUPPORTIVE FOUNDATION LAYERS**
- ✓ Delivers excellent postural comfort
 - ✓ Prevents patient from 'bottoming-out'
 - ✓ High quality, excellent durability

Sonnet Reflex Wedge



The Sonnet Reflex Wedge slightly elevates the head, shoulders and torso to help prevent reflux occurring during sleep. The Sonnet Reflex Wedge can be used while sleeping on side or your without causing any tension in the neck and head.



/ Sonnet Reflex Wedge

CODE	SIZE (MM)
RWS1000	800x880x150/5
RWKS1000	800x1050x150/5
RWEW1000	800x1150x150/5
RWD1000	800x1350x150/5
RWQ1000	800x1530x150/5

Standard Cover:

ADVANTIFLEX

AIROFORM MEMORY IMMERSION INTERFACE

- ✓ Pressure sensitive
- ✓ High temperature stability
- ✓ Visco-elastic cell structure provides excellent immersion and envelopment


SUPPORTIVE FOUNDATION LAYERS

- ✓ Delivers excellent postural comfort
- ✓ Prevents patient from 'bottoming-out'
- ✓ High quality, excellent durability

Sonnet Positioning Wedge



The Sonnet Positioning Wedge is the ideal accessory to achieve 30° side-lying positioning as recommended in clinical best practice publications. Soft and conforming foams provide comfort whilst maintaining sufficient structure and postural support. The unique non-slip base helps maintain position and prevent sliding.



Sonnet Positioning Wedge

CODE	SIZE (MM)
PW1000	280x300x200

Standard Cover:
ADVANTIFLEX

AIROFORM MEMORY IMMERSION INTERFACE

- ✓ Pressure sensitive
- ✓ High temperature stability
- ✓ Visco-elastic cell structure provides excellent immersion and envelopment

SUPPORTIVE FOUNDATION LAYERS

- ✓ Delivers excellent postural comfort
- ✓ Prevents patient from 'bottoming-out'
- ✓ High quality, excellent durability

NON-SLIP BASE MATERIAL

- ✓ Prevents wedge from slipping or moving

Advantiflex Cover

Belgium made polyurethane coated material ensures that Forté Advantiflex Covers deliver consistent care, comfort and durability. Utilising premium quality medical grade fabrics, Forté Advantiflex are a tried and proven medical grade Pressure Relieving Cover.



Premiflex Cover

Forté Premiflex Cover provides the ultimate protection and performance. Superior pressure injury protection is achieved through 4-way stretch non-hammocking and anti-shearing characteristics, Premiflex Covers have a durable double coated polyurethane base that is moisture vapour permeable allowing foam longevity by allowing the core to breath.



Quilted Cotton Cover

Forté Qulited Cotton Cover adds a luxurious touch to any Forte product. Boasting excellent quality and great durability, these hand-made covers facilitate interface immersion and envelopment.





P. 1300 73 11 63

W. fortehealthcare.com.au

E. sales@fortehealthcare.com.au

