




# SONNET PRESSURE OFFLOADING POSITIONING WEDGE



# Sonnet Positioning Wedge



The Sonnet Positioning Wedge is the ideal accessory to achieve 30° side-lying positioning as recommended in clinical best practice publications. Soft and conforming foams provide comfort whilst maintaining sufficient structure and postural support. The unique non-slip base helps maintain position and prevent sliding.



## Sonnet Positioning Wedge

CODE

SIZE (MM)

PW1000

280x300x200

Standard Cover:

**ADVANTIFLEX**

**AIROFORM MEMORY IMMERSION INTERFACE**

- ✓ Pressure sensitive
- ✓ High temperature stability
- ✓ Visco-elastic cell structure provides excellent immersion and envelopment

**SUPPORTIVE FOUNDATION LAYERS**

- ✓ Delivers excellent postural comfort
- ✓ Prevents patient from 'bottoming-out'
- ✓ High quality, excellent durability

**NON-SLIP BASE MATERIAL**

- ✓ Prevents wedge from slipping or moving

# Advantiflex Cover

---

Belgium made polyurethane coated material ensures that Forté Advantiflex Covers deliver consistent care, comfort and durability. Utilising premium quality medical grade fabrics, Forté Advantiflex are a tried and proven medical grade Pressure Relieving Cover.



# Premiflex Cover

---

Forté Premiflex Cover provides the ultimate protection and performance. Superior pressure injury protection is achieved through 4-way stretch non-hammocking and anti-shearing characteristics, Premiflex Covers have a durable double coated polyurethane base that is moisture vapour permeable allowing foam longevity by allowing the core to breath.



# Quilted Cotton Cover

---

Forté Quilted Cotton Cover adds a luxurious touch to any Forte product. Boasting excellent quality and great durability, these hand-made covers facilitate interface immersion and envelopment.





P. 1300 73 11 63

W. [fortehealthcare.com.au](http://fortehealthcare.com.au)

E. [sales@fortehealthcare.com.au](mailto:sales@fortehealthcare.com.au)

