

The logo for FourFront FM4, featuring a stylized icon of three stacked, slanted rectangular blocks to the left of the text "FOURFRONT FM4" in a bold, uppercase, sans-serif font.

FOURFRONT FM4

4-Way Rotational Comfort Mental Health
Mattress

MEDICAL MATTRESS EXPERTS



The Fourfront Range delivers ultimate patient care and support surface versatility within the mental health setting. Fourfront Static Care Mental Health Support Surfaces are the superior choice for demanding mental health ward applications.

Tried and trusted in facilities nationwide, the Fourfront range combines outstanding durability, 4-way rotational ability, a high standard and consistency of care and excellent value.

Fourfront Mattresses incorporate fully sealed and enclosed covers to protect the foam core. All foams used meet fire retardant standards for compliance within Mental Healthcare.



Fourfront FM4 4-Way Rotational Triple Layer Mental Health Mattress



/ Fourfront FM4

CODE	DIMENSIONS (MM)	SIZE
FM4S1000	1980x880x150	Single
FM4SL1000	2030x880x150	Single, Long
FM4KS1000	1980x1050x150	King Single
FM4KSL1000	2030x1050x150	King Single, Long
FM4EW1000	1980x1150x150	Extra-Wide
FM4EWL1000	2030x1150x150	Extra-Wide, Long
FM4D1000	1980x1350x150	Double
FM4DL1000	2030x1350x150	Double, Long
FM4Q1000	1980x1530x150	Queen
FM4QL1000	2030x1530x150	Queen, Long
FM4K1000	1980x1850x150	King
FM4KL1000	2030x1850x150	King, Long
	Customisable	CAN BE MADE TO SUIT ANY BED
THERAPEUTIC WEIGHT RANGE*	150KG	
APPLICATIONS	✔ MENTAL HEALTH WARDS ✔ MENTAL HEALTH FACILITIES ✔ HOME CARE	
WARRANTY**	5-YEAR FOAM CORE, 2-YEAR COVER (12-MONTH WARRANTY IN MENTAL HEALTH APPLICATION)	
STANDARD COVER TYPE	FULLY SEALED & ENCLOSED ADVANTIFLEX	

Features



4-WAY ROTATIONAL TRIPLE LAYERED CORE

- ✓ 4-way rotational core
- ✓ Fire retardant - meets BS5852 source 5 specification
- ✓ Premium grade, high performance foams
- ✓ Profiles & conforms well to Acute Care and Mental Health Beds
- ✓ Foam treated with anti-microbial agent to improve micro-climate and hygiene conditions within the cover



SOFT PATIENT INTERFACE LAYER

- ✓ Enhances patient comfort to deliver superior rest
- ✓ Premium Grade Foam



PREMIUM POSTURAL SUPPORT CENTRE

- ✓ Open cell high quality foam
- ✓ Provides exceptional comfort and postural support for repositioning and patient care
- ✓ Provides exceptional conformity and profiling to medical bed surfaces

SHEAR & FRICTION:



REPOSITIONING:



INFECTION CONTROL:



FIRMNESS:



MOISTURE & MICRO-CLIMATE:



PERFORMANCE & LONGEVITY:




MOBILITY & FUNCTION:



PRESSURE CARE RATING:

LOW-MEDIUM RISK

Standard Custom Options



HINGED STRENGTHENED SIDES

- ✓ Increase lateral stability of support surface
- ✓ Improves patient safety and encourages central patient positioning
- ✓ Make patient transfers easier
- ✓ Unique hinge profile reduces foam fatigue and allow mattress to easily conform to varied bed positions
- ✓ Enables ease of primary patient care



MENTAL HEALTH SPECIFICATION COVER OPTIONS

- ✓ Pressure Care, Vapour Permeable, Polyurethane Material
- ✓ Non-Pressure Care, Damaged Resistant, Waterproof, PVC Material



SIZE

- ✓ Customisable to any size - width, length, height
- ✓ Contact Forté Healthcare

DISCLAIMERS

* Optimal patient weight range for pressure care and comfort. Minimum weight capacities are to be used as a guide only, it may be the situation that patients with a lower weight than the minimum weight threshold will be suited to a surface with a greater minimum weight threshold than that of their own weight for functional and/or individual comfort purposes.

** Warranty refers to any manufacturing or product defects. Due to its nature, foam will soften and its resilience will decrease with usage and age. Contact Forté Healthcare Solutions for further specific information.