




 **FOURFRONT FS4**
Seclusion and Isolation Care Mattress &
Bed Replacement



MEDICAL MATTRESS EXPERTS



The Fourfront Range delivers ultimate patient care and support surface versatility within the mental health setting. Fourfront Static Care Mental Health Support Surfaces are the superior choice for demanding mental health ward applications.

Tried and trusted in facilities nationwide, the Fourfront range combines outstanding durability, 4-way rotational ability, a high standard and consistency of care and excellent value.

Fourfront Mattresses incorporate fully sealed and enclosed covers to protect the foam core. All foams used meet fire retardant standards for compliance within Mental Healthcare.



Fourfront FS4 Seclusion and Isolation Care Mattress/Bed Replacement



/ Fourfront FS

DIMENSIONS (MM)		CODE (FOURFRONT FS4, WITHOUT PILLOW)
1900x900x400	FS41000	
1900x900x500	FS41001	
2000x900x400	FS41002	
2000x900x500	FS41003	
1900x1000x400	FS41004	
1900x1000x500	FS41005	
2000x1000x400	FS41006	
2000x1000x500	FS41007	
THERAPEUTIC WEIGHT RANGE*	N/A	
APPLICATIONS	<input checked="" type="checkbox"/> MENTAL HEALTH WARDS <input checked="" type="checkbox"/> MENTAL HEALTH FACILITIES <input checked="" type="checkbox"/> HOME CARE	
WARRANTY**	5-YEAR FOAM CORE, 4-YEAR COVER (12-MONTH WARRANTY IN MENTAL HEALTH APPLICATION)	
STANDARD COVER TYPE	FULLY SEALED & ENCLOSED HEAVY DUTY PVC	

Features



PATIENT INTERFACE LAYER

- ✓ Enhances patient comfort to deliver superior rest
- ✓ Premium grade, high performance foam
- ✓ Meets BS5852 (fire retardant) source 5 specification

PREMIUM POSTURAL SUPPORT LOWER CORE REGION/FOUNDATION

- ✓ Open cell high quality foam
- ✓ Provides protection against 'bottoming-out' during transfer, repositioning and whilst in the inclined/fowler position
- ✓ Provides exceptional comfort and postural support for repositioning and patient care
- ✓ Meets BS5852 (fire retardant) source 5 specification



OPTIONAL INTEGRATED PILLOW SECTION (FS5)

- ✓ Enhances patient comfort
- ✓ Dual density inverted convolute construction
- ✓ Tamper resistant
- ✓ High performance foam
- ✓ Meets BS5852 source 5 specification

SHEAR & FRICTION:



REPOSITIONING:



INFECTION CONTROL:



FIRMNESS:



MOISTURE & MICRO-CLIMATE:



PERFORMANCE & LONGEVITY:



MOBILITY & FUNCTION:



PRESSURE CARE RATING:

LOW RISK

Standard Custom Options



MENTAL HEALTH SPECIFICATION COVER OPTIONS

- ✓ Pressure Care, Vapour Permeable, Polyurethane Material
- ✓ Non-Pressure Care, Damaged Resistant, Waterproof, PVC Material



SIZE

- ✓ Customisable to any size - width, length, height
- ✓ Contact Forté Healthcare

DISCLAIMERS

* Optimal patient weight range for pressure care and comfort. Minimum weight capacities are to be used as a guide only, it may be the situation that patients with a lower weight than the minimum weight threshold will be suited to a surface with a greater minimum weight threshold than that of their own weight for functional and/or individual comfort purposes.

** Warranty refers to any manufacturing or product defects. Due to its nature, foam will soften and its resilience will decrease with usage and age. Contact Forté Healthcare Solutions for further specific information.