



forté
HEALTHCARE

AIR MATTRESS SURROUNDS

AMS5000 Forté Air Overlay Surrounds

EVERY PRODUCT
CUSTOM MADE

Custom Made to suit any Air Overlay Support Surface

Forté Air Overlay Surrounds provide additional width and/or length to pressure relieving air support surfaces enabling the air support surface to adapt to almost any bed base or configuration. High specification foam enhances the stability of the support surface and provide additional patient and caregiver assistance for positioning and transfers. A firm foam base eliminates the need for a secondary mattress underneath the Air Overlay. Covering materials and construction methods match the clinical performance characteristics of the Air Overlay Support Surface to ensure that pressure relief and patient skin protection is not compromised. Speak with Forté Healthcare or your distributor today to enquire further.

APPLICATIONS

- ✓ Hospital Wards
- ✓ Aged Care Facilities
- ✓ Home Care



OPTIONAL FOAM BASE

- ✓ Provides protection to air cells from sharp protrusions and wear points on the bed base

Air Overlay Surrounds

BRAND	CODE	OVERALL SIZE	
Custom	AMS5002		
Harvest Supreme	AMS5003		
Arjo Nimbus	AMS5004		
Novus Premium 5	AMS5005	King Single	2000x1050
Entrix 5	AMS5006		
Talley Quattro Plus	AMS5007		
Talley Pulsair Choice	AMS5008		
Easy Air 5000	AMS5009		
Custom	AMS5012		
Harvest Supreme	AMS5013		
Arjo Nimbus	AMS5014		
Novus Premium 5	AMS5015	Extra-Wide	2000x1150x200
Entrix 5	AMS5016		
Talley Quattro Plus	AMS5017		
Talley Pulsair Choice	AMS5018		
Easy Air 5000	AMS5019		



PROFILE STABILISING SIDES

- ✓ Enables ease of primary patient care
- ✓ Facilitates patient transfers
- ✓ Unique hinge points allow mattress to easily conform to varied bed positions
- ✓ Encourages central patient positioning

4-WAY STRETCH, VAPOR PERMEABLE, PRESSURE CARE UPPER

- ✓ Breathable
- ✓ Replaces air mattress cover
- ✓ High specification to provide maximum shear relief
- ✓ Available in Premiflex Grade
- ✓ Wipe clean
- ✓ Fire retardant, anti-microbial, resistant to body fluids

HEAVY DUTY PVC BASE

- ✓ Provides enhanced durability
- ✓ Delivers great durability & performance
- ✓ Fire retardant, anti-microbial, resistant to body fluids



COMPATIBLE WITH CPR FUNCTION

- ✓ Custom made to provide access to the CPR function. Note - the access style differs by air mattress brand
- ✓ Cover print ensures visibility for training and in-service purposes



AIR SUPPLY PROTECTION SHEATH

- ✓ Protects against fluid ingress into the cover
- ✓ Easy to clean and manage infection control risks



360° ZIPPER AND WATERFALL INGRESS PROTECTION FLAP

- ✓ Full removable upper cover enables easy installation of air mattress on-site
- ✓ Easy access for internal mattress adjustment/repair
- ✓ Double coated Premiflex Waterfall flap provides protection against fluid ingress. Backing is sealed for ease of cleaning

MEDICAL MATTRESS EXPERTS



P. 1300 73 11 63
W. fortehealthcare.com.au
E. sales@fortehealthcare.com.au