

# forté

HEALTHCARE

Finalé BRP2 Padded Mesh Bed Rail Protectors



# Finalé BRP2 Padded Mesh Bed Rail Protectors



Finalé Padded Mesh Bed Rail Protectors are designed for fitment to full-length bed rails to provide additional patient protection. The integration of a mesh panel provides patient and caregivers with increased visibility and reduces any sense of claustrophobia. Removable and washable, Finalé Padded Mesh Rail protectors are popular with pediatric clients. Custom options include a triple-bar protector and different fabric colours to fit within the care setting.



## ADVANTIFLEX PRESSURE CARE COVER

### ADVANTIFLEX PU1100 POLYURETHANE PRESSURE CARE MATERIAL:

- ✓ Vapour permeable and breathable
- ✓ Anti-shearing
- ✓ Fire retardant, resistant to body fluids and bacteria
- ✓ Wipe-clean
- ✓ Optional: heavy duty PVC material for maximum durability

### MESH INSERT PANEL:

- ✓ Enhances patient and caregiver visibility
- ✓ No compromise of rail functionality and folding
- ✓ Reduces any sense of claustrophobia

### CONSTRUCTION:

- ✓ Closure style: Heavy Duty Zipper fitment - Robust & Secure
- ✓ Robust seams

## 20MM PADDED PROTECTIVE CORE

- ✓ High quality, open cell foam
- ✓ Anti-microbial treatment
- ✓ Padded on both sides of rail for maximum patient protection

## BRP2 Padded Mesh Wrap Around Bed Rail Protectors

COVER	SIZE (MM)	CODE
Advantiflex	1500x400	BRP21000
	1800x400	BRP21001
	2000x400	BRP21002
Heavy Duty PVC	1500x400	BRP25000
	1800x400	BRP25001
	2000x400	BRP25002

## APPLICATIONS

- ✓ Hospitals
- ✓ Aged Care
- ✓ Group Homes
- ✓ Home Care

# MEDICAL MATTRESS EXPERTS



P. 1300 73 11 63  
W. [fortehealthcare.com.au](http://fortehealthcare.com.au)  
E. [sales@fortehealthcare.com.au](mailto:sales@fortehealthcare.com.au)