

forté



HEALTHCARE

Finalé Fall Protection Floor Mat



Finalé Fall Protection Floor Mats

Forté Fall Protection Floor Mats provide impact protection for patients at risk of falling from the bed surface. Care staff can have peace of mind that the Finalé Fall Mat utilises an extremely high impact absorbing foam core for superior patient protection, whilst meeting any infection control and WH&S concerns. Convenient, lightweight, handmade and durable, the Finalé Fall Mats are an essential part of any facility falls-management program.



FULLY SEALED & ENCLOSED HEAVY DUTY PVC2 COVER

HEAVY DUTY, DURABLE PVC1120 MATERIAL:

- ✓ Wipe-clean, Hygienic
- ✓ 100% Waterproof
- ✓ Heavy Duty, Robust & Hardwearing
- ✓ 100% waterproof
- ✓ Anti-Microbial Protection
- ✓ Meets Source 5 Fire Retardant Standardsbody fluids
- ✓ Superior resistance to disinfection chemicals

CONSTRUCTION:

- ✓ Fully sealed and enclosed, no point of entrance to foam core
- ✓ Folds in half for storage
- ✓ Robust seams, Reinforced Critical Stress Points

DUAL DENSITY, HIGH IMPACT ABSORBING FOAM CORE

- ✓ Softens patients impact in event of a fall
- ✓ High quality, durable foam
- ✓ Folds for Storage

FOLDS FOR STORAGE

- ✓ Easy to relocate
- ✓ Compact and space-saving

Finalé Fall Protection Floor Mat - Standard

SIZE (MM)	CODE
1900x750x50	FMS1000

APPLICATIONS

- ✓ Hospitals
- ✓ Aged Care
- ✓ Home Care