



# SONNET PRESSURE CARE CUSHIONS

# Sonnet 2-Core

Pressure Care Cushions



Forté Sonnet 2-Core cushions are manufactured from extremely high performance, medical grade foams and provide sustained patient care and comfort. Graded to deliver both comfort, patient protection and support, the Sonnet cushion range is the ideal choice for facilities seeking low cost, durable and compliant pressure care cushions for low risk patients.



## / Sonnet 2-Core

CODE	SIZE (MM)
CU31100	400x400x70
CU31101	460x400x70
CU31102	460x460x70
CU31103	500x460x70
CU31104	550x460x70
CU31105	550x500x70

Standard Cover:

**ADVANTIFLEX**

**AIROFORM MEMORY MEDIUM IMMERSION  
INTERFACE LAYER**

- ✓ Pressure sensitive
- ✓ High temperature stability
- ✓ Visco-elastic cell structure provides excellent immersion and envelopment

**SUPPORTIVE FOUNDATION LAYER**


- ✓ Delivers excellent postural comfort
- ✓ Prevents patient from 'bottoming-out'
- ✓ High quality, excellent durability

# Sonnet Plush

Triple Core Cushion



Forté Sonnet Plush cushions provide deep profile pressure relief and luxurious patient comfort. 100% made in Forté's New England facility, Sonnet Plush Cushions raise the standard of patient comfort without compromise on clinical effectiveness or suitability for facility usage. Sonnet Plush cushions are also used as a means to provide a higher seating height to assist patient ingress and egress from lounge seating. The Sonnet Plus range is available single or dual sided.



## / Sonnet Plush

CODE	SIZE (MM)
CU51100	400x400x120
CU51101	460x400x120
CU51102	500x460x120
CU51103	500x460x120
CU51104	500x500x120

### Standard Cover:

#### ADVANTIFLEX

#### TRIPLE LAYER DUAL SIDED CORE

- ✓ High quality, open cell foam
- ✓ Anti-microbial treatment
- ✓ Padded on both sides of rail for maximum patient protection
- ✓ 4-way rotational

#### AIROFORM MEMORY EXTRA-SOFT IMMERSION INTERFACE LAYER


- ✓ Pressure sensitive
- ✓ High temperature stability
- ✓ Visco-elastic cell structure provides excellent immersion and envelopment
- ✓ Laminated to both sides of cushion

# Sonnet Bariatric

300kg Triple Core Cushion



The Forté Sonnet Bariatric cushions are manufactured from extremely high performance, medical grade foams and provided sustained patient care and comfort specifically for bariatric patients, up to 300kg's. 100% made in Forté's New England facility, the Sonnet Bariatric cushion has a triple layered foam core graded to deliver both comfort, patient protection and support. The Sonnet cushions are the ideal choice for facilities seeking low cost, durable and compliant pressure care cushions for low risk patients.



## / Sonnet Plush

CODE	SIZE (MM)
CU51100	400x400x120
CU51101	460x400x120
CU51102	500x460x120
CU51103	500x460x120
CU51104	500x500x120

Standard Cover:

**ADVANTIFLEX**

### TRIPLE LAYER DUAL SIDED CORE

- ✓ High quality, open cell foam
- ✓ Anti-microbial treatment
- ✓ Padded on both sides of rail for maximum patient protection
- ✓ 4-way rotational

### AIROFORM MEMORY EXTRA-SOFT IMMERSION INTERFACE LAYER

- ✓ Pressure sensitive
- ✓ High temperature stability
- ✓ Visco-elastic cell structure provides excellent immersion and envelopment
- ✓ Laminated to both sides of cushion

# Advantiflex Cover

---

Belgium made polyurethane coated material ensures that Forté Advantiflex Covers deliver consistent care, comfort and durability. Utilising premium quality medical grade fabrics, Forté Advantiflex are a tried and proven medical grade Pressure Relieving Cover.



# Premiflex Cover

---

Forté Premiflex Cover provides the ultimate protection and performance. Superior pressure injury protection is achieved through 4-way stretch non-hammocking and anti-shearing characteristics, Premiflex Covers have a durable double coated polyurethane base that is moisture vapour permeable allowing foam longevity by allowing the core to breath.



# Quilted Cotton Cover

---

Forté Quilted Cotton Cover adds a luxurious touch to any Forte product. Boasting excellent quality and great durability, these hand-made covers facilitate interface immersion and envelopment.







P. 1300 73 11 63

W. [fortehealthcare.com.au](http://fortehealthcare.com.au)

E. [sales@fortehealthcare.com.au](mailto:sales@fortehealthcare.com.au)

