




SUPER SHEAR INNER COVER



MEDICAL MATTRESS EXPERTS

Super Sheer Anti-Sheering Inner Cover

- Full internal anti-shear lycra stretch inner cover which is an accessory that can be added to any mattress to create a low friction mattress for clients at risk of skin injury due to fragile skin
- Fine stretch lycra based material which is wrinkle free and conforms to varied bed positions
- Provides excellent compensation for shear forces at the patient interface encountered by powered bed articulation, repositioning, and transfer. Instead of the outer cover shearing against the skin during movement, the outer cover shears against the inner cover *INSERT DIAGRAM*
- 100% breathable to maintain microclimatic properties of outer cover and promote inner core breathability
- Excellent bi-directional stretch for uncompromised immersion
- Also provides additional protection to the foam core for increased longevity





P. 1300 73 11 63

W. fortehealthcare.com.au

E. sales@fortehealthcare.com.au