



forté
HEALTHCARE

FINALÉ MEDICAL BED PROTECTORS & ACCESSORIES

BRP3 SLIP-STYLE BED PROTECTORS

Finalé BRP3 Padded Slip-Style Bed Rail Protectors



ADVANTIFLEX PRESSURE CARE COVER

ADVANTIFLEX PU1100 POLYURETHANE PRESSURE CARE MATERIAL:

- ✓ Vapour permeable and breathable
- ✓ Anti-shearing
- ✓ Fire retardant, resistant to body fluids and bacteria
- ✓ Wipe-clean
- ✓ Heavy Duty PVC Inner Material – for maximum durability
- ✓ Optional: heavy duty PVC material for maximum durability

CONSTRUCTION:

- ✓ Closure style: velcro or press - for quick and simple fitment.
- ✓ Robust seams

25MM PADDED PROTECTIVE FOAM CORE

- ✓ High quality, open cell foam
- ✓ Anti-microbial treatment
- ✓ Padded on both sides of rail for maximum patient protection

APPLICATIONS

- ✓ Hospitals
- ✓ Aged Care
- ✓ Group Homes
- ✓ Home Care