



## SONNET 2-CORE PRESSURE CARE CUSHION

# Sonnet 2-Core

Pressure Care Cushions



Forté Sonnet 2-Core cushions are manufactured from extremely high performance, medical grade foams and provide sustained patient care and comfort. Graded to deliver both comfort, patient protection and support, the Sonnet cushion range is the ideal choice for facilities seeking low cost, durable and compliant pressure care cushions for low risk patients.



## / Sonnet 2-Core

CODE	SIZE (MM)
CU31100	400x400x70
CU31101	460x400x70
CU31102	460x460x70
CU31103	500x460x70
CU31104	550x460x70
CU31105	550x500x70

Standard Cover:

**ADVANTIFLEX**

**AIROFORM MEMORY MEDIUM IMMERSION INTERFACE LAYER**

- ✓ Pressure sensitive
- ✓ High temperature stability
- ✓ Visco-elastic cell structure provides excellent immersion and envelopment

**SUPPORTIVE FOUNDATION LAYER**

- ✓ Delivers excellent postural comfort
- ✓ Prevents patient from 'bottoming-out'
- ✓ High quality, excellent durability

# Advantiflex Cover

---

Belgium made polyurethane coated material ensures that Forté Advantiflex Covers deliver consistent care, comfort and durability. Utilising premium quality medical grade fabrics, Forté Advantiflex are a tried and proven medical grade Pressure Relieving Cover.



# Premiflex Cover

---

Forté Premiflex Cover provides the ultimate protection and performance. Superior pressure injury protection is achieved through 4-way stretch non-hammocking and anti-shearing characteristics, Premiflex Covers have a durable double coated polyurethane base that is moisture vapour permeable allowing foam longevity by allowing the core to breath.



# Quilted Cotton Cover

---

Forté Quilted Cotton Cover adds a luxurious touch to any Forte product. Boasting excellent quality and great durability, these hand-made covers facilitate interface immersion and envelopment.





P. 1300 73 11 63

W. [fortehealthcare.com.au](http://fortehealthcare.com.au)

E. [sales@fortehealthcare.com.au](mailto:sales@fortehealthcare.com.au)

