



forté
HEALTHCARE

AIROFORM GEL IMMERSION TECHNOLOGY

AIROFORM IMMERSION TECHNOLOGY

Forté Airoform Immersion technology adapts the latest in Foam expertise with proven Clinical parameters in order to Pressure Injury Prevention and provide optimal Patient Comfort.

Forté Airoform is available in 3 variations, Airoform Memory, Airoform Gel, and Airoform Response.

Each foam technology offers care and comfort unique to patient requirements and can be used individually or blended in order to improve patient care and outcomes.

Airoform Immersion foams are designed to be used at the skin interface level of support surfaces to achieve and improve:

- ✓ Microclimate Conditions
- ✓ Immersion, Pressure Redistribution and Minimising Localised Peak Pressures Damaging to Skin
- ✓ Reduction of Shearing and Related Forces
- ✓ Patient Comfort and Rest
- ✓ Patient Function and Mobility

AIROFORM

gel

Recognised for its patient care attributes as a high performance, cool and conforming foam with Memory Comfort.

MICROCLIMATE

- ✓ Optimises and Moderates microclimatic conditions at the skin interface
- ✓ The most breathable Memory Foam Available
- ✓ Unique Open-Cellular Structure provides superior dissipation of excessive heat
- ✓ Impregnated with Solid-Particle Gel which, when compressed, reduces the reactivity of the foam to applied heat and increases temperature stability

COMFORT

- ✓ Available in Ultra Soft

SKIN INTEGRITY AND PRESSURE RELIEF

- ✓ Provides Excellent Pressure Redistribution and extremely low interface pressure
- ✓ Maximises the Contact Between the Patient and the Surface, Decreasing Interface Pressure
- ✓ Offers excellent heat stability, provides sustained immersion reactivity of the foam to applied heat and increases temperature stability

CONTRAINDICATIONS

- ✓ Must be used in conjunction with Forté Anti-Shearing Inner Cover Technology
- ✓ May not suit 'cold sleeper'



| | |
|-----------------------------|--|
| Shear & Friction | 0 Respective to application 5 |
| Repositioning | 0 Respective to application 5 |
| Infection Control | Respective to application |
| Firmness | 0 Respective to application 5 |

| | |
|------------------------------------|--|
| Moisture & Microclimate | 0 Respective to application 5 |
| Performance & Longevity | 0 Respective to application 5 |
| Mobility & Function | 0 Respective to application 5 |
| Pressure Care Rating | Respective to application |

| Code | Dimensions (L x W x D mm) | Size |
|---|---------------------------------------|--|
| AIROFORM GEL IMMERSION INTERFACE LAYER, MATTRESS CUSTOMISATION | | |
| GEL | Respective to Selected Mattress | CUSTOMISABLE |
| AIROFORM GEL IMMERSION INTERFACE LAYER W/ANTI-SHEARING INNER COVER, MATTRESS CUSTOMISATION | | |
| GEL-IC | Respective to Selected Mattress | CUSTOMISABLE |
| AIROFORM GEL IMMERSION OVERLAY – FOAM & STOCKINETTE COVER | | |
| OLGELS1000 | 1980x880x50 | Single |
| OLGELSL1000 | 2030x880x50 | Single, Long |
| OLGELKS1000 | 1980x1050x50 | King Single |
| OLGELKSL1000 | 2030x1050x50 | King Single, Long |
| OLGELEW1000 | 1980x1150x50 | Extra-Wide |
| OLGELEWL1000 | 2030x1150x50 | Extra-Wide, Long |
| OLGELD1000 | 1980x1350x50 | Double |
| OLGELDL1000 | 2030x1350x50 | Double, Long |
| OLGELQ1000 | 1980x1530x50 | Queen |
| OLGELQL1000 | 2030x1530x50 | Queen, Long |
| OLGELK1000 | 1980x1850x50 | King |
| OLGELKL1000 | 2030x1850x50 | King, Long |
| | CUSTOMISABLE Available with Covers | To suit any bed platform / mattress |

Therapeutic Weight Range • Respective to its Application

Applications

- Hospital Wards
- Aged Care
- Palliative Care
- Homecare

Warranty

- 5 Years
- Respective to its Application

Standard Cover • Respective to its Application



GLOBAL TECHNOLOGY
AUSTRALIAN MADE

All mattresses are made by
Forté Healthcare in Armidale NSW
in accordance with all Australian
manufacturing
and clinical standards.

