




forté
HEALTHCARE

SONNET PRESSURE OFFLOADING HEEL DEVICES

Sonnet Heel Wedge



The Sonnet Heel Wedge is used for leg elevation and heel suspension for effective pressure relief to critical heel areas. Manufactured using premium immersion foams and a pressure care covering, the Forté Heel Wedge is suited for a wide range of patients and ideal for performance and longevity.



/ Sonnet Heel Wedge

| CODE | SIZE (MM) |
|--------|------------|
| HW1000 | 280mm wide |
| HW1001 | 400mm wide |
| HW1002 | 680mm wide |

Standard Cover:

ADVANTIFLEX

AIROFORM MEMORY IMMERSION INTERFACE

- ✓ Pressure sensitive
- ✓ High temperature stability
- ✓ Visco-elastic cell structure provides excellent immersion and envelopment

SUPPORTIVE FOUNDATION LAYERS

- ✓ Delivers excellent postural comfort
- ✓ Prevents patient from 'bottoming-out'
- ✓ High quality, excellent durability

Advantiflex Cover

Belgium made polyurethane coated material ensures that Forté Advantiflex Covers deliver consistent care, comfort and durability. Utilising premium quality medical grade fabrics, Forté Advantiflex are a tried and proven medical grade Pressure Relieving Cover.



Premiflex Cover

Forté Premiflex Cover provides the ultimate protection and performance. Superior pressure injury protection is achieved through 4-way stretch non-hammocking and anti-shearing characteristics, Premiflex Covers have a durable double coated polyurethane base that is moisture vapour permeable allowing foam longevity by allowing the core to breath.



Quilted Cotton Cover

Forté Quilted Cotton Cover adds a luxurious touch to any Forte product. Boasting excellent quality and great durability, these hand-made covers facilitate interface immersion and envelopment.





GLOBAL TECHNOLOGY
AUSTRALIAN MADE

All mattresses are made by
Forté Healthcare in Armidale NSW
in accordance with all Australian
manufacturing
and clinical standards.

