



MEDICAL MATTRESS EXPERTS



The Fourfront Range delivers ultimate patient care and support surface versatility within the mental health setting. Fourfront Static Care Mental Health Support Surfaces are the superior choice for demanding mental health ward applications.

Tried and trusted in facilities nationwide, the Fourfront range combines outstanding durability, 4-way rotational ability, a high standard and consistency of care and excellent value.

Fourfront Mattresses incorporate fully sealed and enclosed covers to protect the foam core. All foams used meet fire retardant standards for compliance within Mental Healthcare.



Fourfront FM6 4-Way Rotational Ventraflow Mental Health Mattress



SHEAR & FRICTION:



REPOSITIONING:



INFECTION CONTROL:



FIRMNESS:



MOISTURE & MICRO-CLIMATE:



PERFORMANCE & LONGEVITY:



MOBILITY & FUNCTION:



PRESSURE CARE RATING:

LOW-MEDIUM RISK

/ Fourfront FM6

CODE	DIMENSIONS (MM)	SIZE
FM6S1000	1980x880x150	Single
FM6SL1000	2030x880x150	Single, Long
FM6KS1000	1980x1050x150	King Single
FM6KSL1000	2030x1050x150	King Single, Long
FM6EW1000	1980x1150x150	Extra-Wide
FM6EWL1000	2030x1150x150	Extra-Wide, Long
FM6D1000	1980x1350x150	Double
FM6DL1000	2030x1350x150	Double, Long
FM6Q1000	1980x1530x150	Queen
FM6QL1000	2030x1530x150	Queen, Long
FM6K1000	1980x1850x150	King
FM6KL1000	2030x1850x150	King, Long
	Customisable	CAN BE MADE TO SUIT ANY BED
THERAPEUTIC WEIGHT RANGE*	150KG	
APPLICATIONS	✓ MENTAL HEALTH WARDS ✓ MENTAL HEALTH FACILITIES ✓ HOME CARE	
WARRANTY**	5-YEAR FOAM CORE, 4-YEAR COVER (12-MONTH WARRANTY IN MENTAL HEALTH APPLICATION)	
STANDARD COVER TYPE	FULLY SEALED & ENCLOSED PREMIFLEX	

Features



UNIQUE, 5 ZONED CORE

- ✓ High specification foam for longevity and performance
- ✓ Provides graded postural support and patient comfort
- ✓ Anti-microbial treatment applied to foam maximizes longevity and core micro-climate
- ✓ Conforms to EPUAP, NPUAP & PPPIA Guidelines for quality and comfort
- ✓ Fire retardant - meets BS5852 source 5 specification
- ✓ Superior ability to profile and conform to medical bed surfaces



AIROFORM RESPONSE IMMERSION INTERFACE LAYER

- ✓ A particularly responsive Ultra-High Performance elastic foam designed to optimise patient mobility, function and enhance repositioning
- ✓ Open-cell structure provides low-shear surface, good pressure redistribution, effective immersion, enhanced breathability and superior heat dissipation
- ✓ Alleviates localised pressure points
- ✓ Enhances contact with patient for distribution of interface pressure
- ✓ High temperature stability to deliver consistent immersion, performance and comfort
- ✓ Provides exceptional comfort
- ✓ Greatly reduces risk of impression marks on patient skin
- ✓ Firmness: Medium



VENTRAFLOW® INNER CORE TECHNOLOGY

- ✓ Profile support design conforms to patient contact and immerses critical risk pressure areas i.e. bony prominences
- ✓ Delivers postural support for patient reposition
- ✓ Internal air circulation provides over 3 times greater internal air flows minimizing internal heat entrapment and maximum moisture evacuation
- ✓ Triple support region to deliver enhanced immersion and support for critical body regions i.e. heel/scapular/sacrum
- ✓ Enhances core stability by providing additional support where it is needed



PREMIUM POSTURAL SUPPORT MID-CORE REGION

- ✓ Graded support and patient comfort in supine, prone, lateral and fowler positions
- ✓ Provides protection against 'bottoming-out' during transfer, repositioning and whilst in the inclined/fowler position
- ✓ Provides exceptional comfort, support for repositioning and patient care

HINGED STRENGTHENED SIDES

- ✓ Increase lateral stability of support surface
- ✓ Improves patient safety and encourages central patient positioning
- ✓ Make patient transfers easier
- ✓ Unique hinge profile reduces foam fatigue and allow mattress to easily conform to varied bed positions
- ✓ Enables ease of primary patient care

Standard Custom Options



HINGED STRENGTHENED SIDES

- ✓ Increase lateral stability of support surface
- ✓ Improves patient safety and encourages central patient positioning
- ✓ Make patient transfers easier
- ✓ Unique hinge profile reduces foam fatigue and allow mattress to easily conform to varied bed positions
- ✓ Enables ease of primary patient care



MENTAL HEALTH SPECIFICATION COVER OPTIONS

- ✓ Pressure Care, Vapour Permeable, Polyurethane Material
- ✓ Non-Pressure Care, Damaged Resistant, Waterproof, PVC Material



SIZE

- ✓ Customisable to any size - width, length, height
- ✓ Contact Forté Healthcare

DISCLAIMERS

* Optimal patient weight range for pressure care and comfort. Minimum weight capacities are to be used as a guide only, it may be the situation that patients with a lower weight than the minimum weight threshold will be suited to a surface with a greater minimum weight threshold than that of their own weight for functional and/or individual comfort purposes.

** Warranty refers to any manufacturing or product defects. Due to its nature, foam will soften and its resilience will decrease with usage and age. Contact Forté Healthcare Solutions for further specific information.



GLOBAL TECHNOLOGY
AUSTRALIAN MADE

All mattresses are made by
Forté Healthcare in Armidale NSW
in accordance with all Australian
manufacturing
and clinical standards.

