



forté
HEALTHCARE

FOURFRONT FS5

Seclusion and Isolation Care All-in-One Mattress



The Fourfront Range delivers ultimate patient care and support surface versatility within the mental health setting. Fourfront Static Care Mental Health Support Surfaces are the superior choice for demanding mental health ward applications.

Tried and trusted in facilities nationwide, the Fourfront range combines outstanding durability, 4-way rotational ability, a high standard and consistency of care and excellent value.

Fourfront Mattresses incorporate fully sealed and enclosed covers to protect the foam core. All foams used meet fire retardant standards for compliance within Mental Healthcare.



Fourfront FS5 Seclusion and Isolation Care Mattress/Bed Replacement with integrated Pillow



/ Fourfront FS

DIMENSIONS (MM)	CODE FOURFRONT FS5, (WITH PILLOW)
1900x900x400	FS51000
1900x900x500	FS51001
2000x900x400	FS51002
2000x900x500	FS51003
1900x1000x400	FS51004
1900x1000x500	FS51005
2000x1000x400	FS51006
2000x1000x500	FS51007
THERAPEUTIC WEIGHT RANGE*	N/A
APPLICATIONS	✔ MENTAL HEALTH WARDS ✔ MENTAL HEALTH
FACILITIES	✔ HOME CARE
WARRANTY** WARRANTY IN MENTAL HEALTH APPLICATION)	5-YEAR FOAM CORE, 4-YEAR COVER (12-MONTH
STANDARD COVER TYPE	FULLY SEALED & ENCLOSED HEAVY DUTY PVC

Features



PATIENT INTERFACE LAYER

- ✓ Enhances patient comfort to deliver superior rest
- ✓ Premium grade, high performance foam
- ✓ Meets BS5852 (fire retardant) source 5 specification

PREMIUM POSTURAL SUPPORT LOWER CORE REGION/FOUNDATION

- ✓ Open cell high quality foam
- ✓ Provides protection against 'bottoming-out' during transfer, repositioning and whilst in the inclined/fowler position
- ✓ Provides exceptional comfort and postural support for repositioning and patient care
- ✓ Meets BS5852 (fire retardant) source 5 specification



INTEGRATED PILLOW SECTION

- ✓ Enhances patient comfort
- ✓ Dual density inverted convolute construction
- ✓ Tamper resistant
- ✓ High performance foam
- ✓ Meets BS5852 source 5 specification

SHEAR & FRICTION:



REPOSITIONING:



INFECTION CONTROL:



FIRMNESS:



MOISTURE & MICRO-CLIMATE:



PERFORMANCE & LONGEVITY:



MOBILITY & FUNCTION:



PRESSURE CARE RATING:


LOW RISK

Standard Custom Options



MENTAL HEALTH SPECIFICATION COVER OPTIONS

- ✓ Pressure Care, Vapour Permeable, Polyurethane Material
- ✓ Non-Pressure Care, Damaged Resistant, Waterproof, PVC Material



SIZE

- ✓ Customisable to any size - width, length, height
- ✓ Contact Forté Healthcare

DISCLAIMERS

* Optimal patient weight range for pressure care and comfort. Minimum weight capacities are to be used as a guide only, it may be the situation that patients with a lower weight than the minimum weight threshold will be suited to a surface with a greater minimum weight threshold than that of their own weight for functional and/or individual comfort purposes.

** Warranty refers to any manufacturing or product defects. Due to its nature, foam will soften and its resilience will decrease with usage and age. Contact Forté Healthcare Solutions for further specific information.



GLOBAL TECHNOLOGY
AUSTRALIAN MADE

All mattresses are made by
Forté Healthcare in Armidale NSW
in accordance with all Australian
manufacturing
and clinical standards.

