



forté
HEALTHCARE

AIROFORM RESPONSE IMMERSION TECHNOLOGY

AIROFORM IMMERSION TECHNOLOGY

Forté Airoform Immersion technology adapts the latest in Foam expertise with proven Clinical parameters in order to Pressure Injury Prevention and provide optimal Patient Comfort.

Forté Airoform is available in 3 variations, Airoform Memory, Airoform Gel, and Airoform Response.

Each foam technology offers care and comfort unique to patient requirements and can be used individually or blended in order to improve patient care and outcomes.

Airoform Immersion foams are designed to be used at the skin interface level of support surfaces to achieve and improve:

- ✓ Microclimate Conditions
- ✓ Immersion, Pressure Redistribution and Minimising Localised Peak Pressures Damaging to Skin
- ✓ Reduction of Shearing and Related Forces
- ✓ Patient Comfort and Rest
- ✓ Patient Function and Mobility

AIROFORM RESPONSE

Ultra-High Performance open-cell foam designed to optimise patient mobility, function and enhance repositioning

MICROCLIMATE

- ✓ Open-Cellular Structure for superior heat dissipation

COMFORT


- ✓ Available in Soft and Medium
- ✓ Provides Excellent Postural Support

SKIN INTEGRITY AND PRESSURE RELIEF

- ✓ Supports repositioning and bed mobility through its high responsiveness
- ✓ High Elasticity providing good pressure redistribution
- ✓ Maximises the Contact Between the Patient and the Surface, Decreasing Interface Pressure
- ✓ Offers excellent heat stability, provides sustained immersion

CONTRAINDICATIONS

- ✓ Airoform Response Medium - When used as interface layer together with high risk surface (i.e. Sovereign Mattress), thickness should not exceed 15-20mm
- ✓ May not suit 'cold sleeper'



Shear & Friction	0 Respective to application 5
Repositioning	0 Respective to application 5
Infection Control	Respective to application
Firmness	Available in Soft & Medium

Moisture & Microclimate	0 Respective to application 5
Performance & Longevity	0 Respective to application 5
Mobility & Function	0 Respective to application 5
Pressure Care Rating	Respective to application

Code**Dimensions** (L x W x D mm)**Size****AIROFORM RESPONSE SOFT IMMERSION INTERFACE LAYER, MATTRESS CUSTOMISATION**

ARS	Respective to Selected Mattress	CUSTOMISABLE
-----	---------------------------------	--------------

AIROFORM RESPONSE MEDIUM IMMERSION INTERFACE LAYER, MATTRESS CUSTOMISATION

ARM	Respective to Selected Mattress	CUSTOMISABLE
-----	---------------------------------	--------------

AIROFORM RESPONSE SOFT IMMERSION OVERLAY – FOAM ONLY

OLARSS1000	1980x880x50	Single
OLARSSL1000	2030x880x50	Single, Long
OLARSKS1000	1980x1050x50	King Single
OLARSKSL1000	2030x1050x50	King Single, Long
OLARSEW1000	1980x1150x50	Extra-Wide
OLARSEWL1000	2030x1150x50	Extra-Wide, Long
OLARSD1000	1980x1350x50	Double
OLARSDL1000	2030x1350x50	Double, Long
OLARSQ1000	1980x1530x50	Queen
OLARSQL1000	2030x1530x50	Queen, Long
OLARSK1000	1980x1850x50	King
OLARSKL1000	2030x1850x50	King, Long
	CUSTOMISABLE Available with Covers	To suit any bed platform / mattress

AIROFORM RESPONSE MEDIUM IMMERSION OVERLAY – FOAM ONLY

OLARMS1000	1980x880x50	Single
OLARMSL1000	2030x880x50	Single, Long
OLARMKS1000	1980x1050x50	King Single
OLARMKSL1000	2030x1050x50	King Single, Long
OLARMEW1000	1980x1150x50	Extra-Wide
OLARMEWL1000	2030x1150x50	Extra-Wide, Long
OLARMD1000	1980x1350x50	Double
OLARMDL1000	2030x1350x50	Double, Long
OLARMQ1000	1980x1530x50	Queen
OLARMQL1000	2030x1530x50	Queen, Long
OLARMK1000	1980x1850x50	King
OLARMKL1000	2030x1850x50	King, Long
	CUSTOMISABLE Available with Covers	To suit any bed platform / mattress

Therapeutic Weight Range

- Respective to its Application
-

Applications

- Hospital Wards
 - Aged Care
 - Palliative Care
 - Homecare
-

Warranty

- 5 Years
- Respective to its Application
-

Standard Cover

- Respective to its Application
-



GLOBAL TECHNOLOGY
AUSTRALIAN MADE

All mattresses are made by
Forté Healthcare in Armidale NSW
in accordance with all Australian
manufacturing
and clinical standards.

